



neiform MED rapid



Rapid disinfectant for medical devices

Ready-for-use solution

Fields of application:

- Specific and rapid disinfection of surfaces of medical devices such as surgical tables in hospitals as well as in doctors' and dentists' practices

Performance spectrum:

- Bactericidal, yeasticidal, fungicidal and tuberculocidal activity confirmed in accordance with VAH¹ and EN methods²; activity against enveloped viruses confirmed according to the test guidelines of the RKI/DVV³ by certification
- Also limited spectrum virucidal activity (active against noro and adenovirus) and active against rotavirus
- The application in the bucket wipes system Kimtech* Wettask* DS of Kimberly-Clark fulfils the requirements of the current VAH recommendation⁴
- The surfaces to be disinfected must be compatible with alcoholic disinfectants
- Not suitable for disinfecting acrylic glass (perspex)
- VAH-listed
- Included in the IHO⁵ list of virucidal products

Special properties:

- Wide activity spectrum, even with short contact times
- Rapid and stain-free drying, no rinsing with fresh water necessary
- Free of aldehydes, colourants and perfumes

Application and dosage:

Apply neiform MED rapid undiluted on clean and dry surfaces or on a clean disposable cloth and wipe surfaces. Surfaces must be completely wetted with neiform MED rapid. To achieve the desired activity, the below-mentioned contact times must be strictly adhered to.



neiform MED rapid can also be used with the refillable wiping system Kimtech* Wettask* DS of Kimberly-Clark. Place the roll of wipes into the wipes dispenser system. Fill 2.25 l neiform MED rapid undiluted into the bucket or 0.75 l into the box. Tightly close the wipes dispenser system. After 60 minutes the wipes are completely soaked and ready for use. Simply pull out the soaked wipes through the opening of the dispenser in the required amount. After removal tightly close the dispenser system.

The soaked wipes can be used for up to 28 days when the wiping system is kept tightly closed.

The VAH recommends thorough disinfecting cleaning and drying after depletion and before replenishing the dispenser system. The wipes system Kimtech* Wettask* DS of Kimberly-Clark can be reprocessed manually as well as thermally (thermostable up to 100 °C) or chemothermally in a washer disinfectant.

	Clean conditions, 20 °C (cleaned surfaces)
bactericidal (VAH, EN 13727, EN 16615)	undiluted, 1 min
yeasticidal (VAH, EN 13624, EN 16615)	undiluted, 1 min
fungicidal (EN 13624)	undiluted, 20 min
tuberculocidal (EN 14348)	undiluted, 30 sec
Active against enveloped viruses (incl. HIV, HBV, HCV) (EN 14476, RKI/DVV)	undiluted, 15 sec
Limited virucidal activity (adenovirus, norovirus, incl. all enveloped viruses) (EN 14476)	undiluted, 2 min
Active against norovirus (MNV ⁶) (EN 14476)	undiluted, 2 min
Active against adenovirus (EN 14476)	undiluted, 1 min
Active against rotavirus (EN 14476)	undiluted, 1 min
VAH list	undiluted, 5 min



neoform MED rapid

Due to its alcohol content, neoform MED rapid must not be applied in too large amounts (apply max. 50 ml per m² on the surfaces to be treated and max. 100 ml per m² base area). Please thoroughly read and observe the hazard and precautionary statements.

General instructions for use:

- For professional use only.
- For economical and controlled dosage the use of manual dosing aids is recommended. Please contact us.
- It is generally recommended to wear gloves when working with disinfectants.
- Reprocessing should comply with all ordinances pursuant to the Medical Device Directive and should be performed with appropriate validated processes.
- Please observe the medical device manufacturer's recommendations for use.
- Do not mix with other products.

Expert reports:

The disinfecting activity has been confirmed by certification. We will be pleased to provide certificates on request.

Technical data:

Viscosity	< 10 mPa s (undiluted, 20 °C)
Density	0.9 g/cm ³ (20 °C)
Flash point	27 °C according to DIN 22791

Ingredients:


Active substances in 100 g:
35.0 g propan-1-ol,
25.0 g ethanol

CE-mark:  0297

neoform MED rapid conforms to the European directive 93/42/EC, Annex I concerning medical devices.

Storage information:

Always store at a temperature between -15 °C and 25 °C. Keep container tightly closed.

Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see EC safety data sheets. These are available at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

MB 4070/3-4 Date of Issue: 08/2017

- 1 Verbund für Angewandte Hygiene [Association of Applied Hygiene, Germany]
- 2 European test methods
- 3 according to the guidelines of the Robert Koch-Institute (RKI) and the Deutschen Vereinigung zur Bekämpfung von Viruskrankheiten (DVV) [German Association for the Control of Virus Diseases]
- 4 Empfehlung zur Kontrolle kritischer Punkte bei der Anwendung von Tuchspendensystemen im Vortränksystem für die Flächendesinfektion (Hyg Med 2012; 37-11) [Recommendation for checking critical points when using wipes dispenser systems in a pre-soak system for surface disinfection]
- 5 Industrieverband Hygiene und Oberflächenschutz [Industry Association Surface Protection and Hygiene, Germany]
- 6 Murine norovirus



With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.